

Gamma Rays

OBJECT:

- To understand the various interactions of gamma rays with matter.
- To calibrate a gamma ray scintillation spectrometer, using gamma rays of known energy, and use it to measure the energy of an “unknown” gamma ray.
- To use positron annihilation radiation to determine the mass of the electron and to observe correlated gamma rays.

READINGS: The lab manual “Experiments in Nuclear Science” AN34, 1984 EG&G ORTEC provides an excellent hands-on discussion of the background and techniques for a number of undergraduate level nuclear experiments. The equipment described resembles, with some variation, the equipment available in the laboratory. Additional readings are given at the end of this write-up.

APPARATUS: NaI:Tl scintillator and photomultiplier tube detector with integrated preamplifier (2), high voltage power supply, Canberra model 2000 power supply, NIM bin, Canberra 2015A amplifier/single channel analyzer module (2), Rutgers P1075 scaler/timer, Canberra model 1446 coincidence module, multichannel analyzer PC board in computer, Maestro Analyzer software, monitor.

BACKGROUND: In this experiment you will study the radioactive decay of a nucleus by detecting gamma rays emitted consequent to the decay. Gamma ray detection is a slightly complicated, multi-step process: the gamma ray enters a NaI:Tl scintillator crystal where it produces a rapidly moving free electron that, in turn, loses its energy by excitation of the ions in its path as it travels through the crystal. This excitation energy is given off in various ways, one being emission of visible light (fluorescence). Thus a single high energy gamma ray entering the scintillator produces a flash of low energy photons.

These photons are directed to the photosensitive surface of a photomultiplier tube, where they eject electrons via the photoelectric effect. The electrons are collected in the photomultiplier and amplified to yield a current pulse, which is converted to a voltage pulse whose height is proportional to the number of photoelectrons and is thus proportional to the number of photons reaching the tube, which in turn is proportional to the initial energy of the fast electron.

When a radioactive source is placed near the scintillator, the photomultiplier produces a series of pulses, each corresponding to the decay of a single nucleus. The amplitude of each pulse is related to the energy of the electron freed by the gamma ray. These pulses are studied using either a single- or multi-channel

analyzer. A single channel analyzer (SCA) counts the number of voltage pulses whose height falls within a given (adjustable) window of values, while a multi-channel analyzer (MCA) sorts the pulses according to height and counts the number in each window to give a spectral (energy) distribution of the fast electrons. Figure 1 shows a typical MCA spectrum. In order to relate this spectrum to the nuclear decay, we need to understand how gamma rays interact with matter.




Figure 1. Na(Tl) Spectrum for Cs^{137}

Interaction mechanisms for gamma rays with matter: When entering a crystal, gamma rays produce fast electrons by three different processes: the photoelectric effect, the Compton effect (Compton scattering) and pair production. It is these fast electrons, which give rise to scintillations, not the gamma ray. The observed spectral distribution will thus depend on the detailed interaction process of the gamma rays in the crystal.

Consider a beam of mono-energetic gamma rays striking the scintillator. For our purposes the most important energy loss mechanism is the **photoelectric effect**. When a gamma ray strikes an ion in the crystal, it is absorbed and all of its energy is transferred to one of the bound electrons, which is freed and moves rapidly through the crystal. Since the energy of the gamma ray (typically about 0.5 MeV) is much greater than the binding energy of the electron of the ion (typically 10 to 100 eV), the energy of the freed electron can be considered equal to that of the incoming gamma. (Especially since the energy resolution of the detector is only about 10%.) Thus the photoelectric effect results in a peak, called the **photopeak**, in the photomultiplier spectrum at an energy equal to that of the

incoming gamma ray.

In **Compton scattering**, the gamma ray is not absorbed, but rather scattered through an angle θ by an electron, which recoils and carries away some of the gamma ray's energy E . (The scattered gamma ray then escapes from the scintillator; the probability that a gamma ray Compton scatters in a typical size scintillator is quite small (1% to 10%), which means you are unlikely to detect a gamma ray that has undergone two Compton scatterings.) The gamma ray's initial wavelength is $\lambda = hc/E = 1240/E$ nm, where E is in eV. The change in wavelength is:

$$\Delta\lambda = h/mc (1 - \cos \theta) = 0.00243 (1 - \cos \theta) \text{ nm}, \quad (1)$$

where h is Planck's constant, m is the mass of the electron and c is the speed of light. From this equation you can see that the energy of the scattered electron, which is the energy loss of the gamma ray, will vary from zero (when $\theta = 0^\circ$) to a maximum corresponding to a wavelength shift of 0.00486 nm (when $\theta = 180^\circ$). This maximum energy is called the **Compton edge**. The energy distribution of Compton scattered electrons is essentially a constant. So the Compton spectrum produced by a photomultiplier tube is an almost flat plateau from zero energy up to the Compton edge where it drops off sharply (at a rate limited by the energy resolution of the tube).

Consider, for example, a 622 keV gamma ray from a Cs^{137} decay. Its initial wavelength is 0.00199 nm. If it is Compton scattered through $\theta = 180^\circ$ the wavelength becomes $0.00199 + 0.00486 = 0.00685$ nm corresponding to an energy of 181 keV. The lost energy, $622 - 181 = 441$ keV, is transferred to the electron from which it scattered and is observed as a Compton edge at 441 keV. For smaller scattering angles, the transferred energy will be less, extending down to zero for $\theta = 0^\circ$.

The discussion above refers to gamma rays that are Compton scattered by electrons within the scintillator. It is also possible for a gamma ray to be Compton scattered into the scintillator from an interaction outside the scintillator. In this case the observed signal is from the scattered gamma and not from the recoiling electron. The scattered gamma ray could then be detected through the photoelectric effect. For Cs^{137} the Compton scattered gamma rays will have energies ranging from 181 keV up to the full 622 keV. However, because of the geometry of the detector, most of the gamma rays scattered into the scintillator will have been scattered through a large value of θ . But $\cos \theta$ varies only slowly with θ for θ near 180° , which means [see Eq. (1)] that these gamma rays will all have energies near 181 keV. The resulting energy peak is called the **backscatter peak**. It can be enhanced by placing a sheet of lead around

the outside of the scintillator.

The third interaction mechanism is **pair production**. If the incoming gamma ray energy is above $1.02 \text{ MeV} = 2mc^2$, the rest mass of an electron-positron pair, the gamma ray can spontaneously create an electron-positron pair and be totally absorbed. If both the electron and positron lose all of their kinetic energy while still in the scintillator, they would produce a photomultiplier pulse corresponding to an energy $2mc^2$ below the gamma ray energy E . (Of course, either might escape the crystal after partial loss of kinetic energy.)

But the spectrum is actually more complicated since if the positron has been slowed down and stopped in the crystal, it will annihilate with an electron, emitting two gamma rays, each of energy mc^2 . One, or both, of these gamma rays may be absorbed in the crystal, and thus contribute to the height of the photomultiplier output pulse. Pair production therefore produces a “**full energy**” peak (E), a “**one-escape peak**” ($E - mc^2$) and “**two escape peak**” ($E - 2mc^2$), depending upon whether both annihilation photons are absorbed in the scintillator or one, or both, escapes.

The final question to consider is that of the relative importance of the three interaction mechanisms, which depend in different ways upon the energy of the gamma ray. For low energy rays, the photoelectric effect predominates. Since the photopeak directly yields the energy of the gamma ray, most scintillators are designed to maximize the photopeak. In the NaI:Tl scintillator you will use, a small amount of the heavy metal thallium is added for this purpose when the crystal is grown. (The strength of the photoelectric effect depends strongly on the number of electrons bound to the ion.) As E increases, the photoelectric absorption decreases rapidly, while the Compton scattering decreases much more slowly and predominates above several hundred keV. The absorption coefficient for pair production rises rapidly above the threshold $E = 1.02 \text{ MeV}$ and exceeds the Compton scattering, while photoelectric absorption becomes negligible.

Energy resolution: A NaI:Tl detector has an energy resolution of only about 10%. [More advanced detectors such as high purity germanium detectors have as much as 30 times greater resolution, but are much more expensive.] When a beam of mono-energetic gamma rays strikes the scintillator, there is a fluctuation from gamma ray to gamma ray in the height of the voltage pulse from the photomultiplier, which shows up as a broadening of the photopeak. The pulse height variation is chiefly due to statistical fluctuations in the number of electrons emitted at the cathode of the photomultiplier when a flash of photons arrives from the scintillator, but is also due to the occasional escape of electrons, X-rays or gamma rays from the crystal, all of which depend on how large the NaI crystal is. The **full width of a peak at half its maximum height** (FWHM) is a

convenient measure of the **resolution** of the instrument.

Nuclear decay: You will study the decay of several light nuclei that proceed by beta decay -- the emission of an electron or positron (and the associated neutrino). These decays, Figure 2, leave the daughter nucleus in an excited state that then returns to the ground state by the emission of one or more gamma rays. The range of the beta particles is so short that they cannot penetrate into the




Figure 2. Principal Decay Modes of Cesium 137, Cobalt 60, and Sodium 22.

scintillator. So you will only be detecting the gamma rays. Table 1 gives the energies of the gamma rays emitted by the daughter nuclei for the three nuclei you will be studying. The Ba K X-ray occurs about 4% of the time when instead of emitting a gamma ray, the Ba nucleus de-excites by transferring its energy to one of the innermost electrons (K shell) that overlaps the nucleus. This gives the electron enough energy to escape from the atom. One of the outer electrons then drops into the hole left by the escaping electron and emits its excess energy as an X-ray. Although the escaping electron cannot penetrate into the scintillator, the

K X-ray does and is efficiently detected.

Parent nucleus	Lifetime (years)	Decay mode	Daughter nucleus	Gamma energy (MeV)
Na ²²	2.605	e ⁺	Ne ²²	1.2746
Co ⁶⁰	5.272	e ⁻	Ni ⁶⁰	1.1732
				1.3325
Cs ¹³⁷	30.17	e ⁻	Ba ¹³⁷	0.6616
				0.0322*

*Ba k X-ray

Table 1. Nuclear Decay Data
(from: Handbook of Chemistry and Physics)

γ - γ Angular Correlations: When the nuclei in a radioactive sample decay the individual decays are unrelated and the emitted particles (gamma rays) are emitted isotropically in all directions. We will study one case where the gamma rays from a sample are correlated -- the annihilation radiation of positrons. When a Na²² nucleus emits a positron, the positron rapidly slows and then annihilates, upon encountering an electron. Conservation of momentum requires that the two gamma rays that are emitted travel in exactly opposite directions. You will use two scintillation counters and a coincidence module to study this angular correlation.




Figure 3: Block Diagram for Scintillation Counter Apparatus

PROCEDURE: Warning. *The experiment involves high voltages; do not turn on any switches or knobs until you have become familiar with the apparatus. The high voltage applied to the phototube should never exceed 1200 volts. Observe caution at all times with the scintillation detector since the crystal and photomultiplier are fragile and expensive.*

The apparatus is for the most part pre-wired and very simple to operate. The block diagram is shown in Fig. 3. The scintillation counter consists of a NaI:TI crystal glued to a photomultiplier (PM) tube. The PM tube is powered by an adjustable high voltage supply (~1000 V). The supply is a vacuum tube model and should be set to “standby” for about 1 minute before turning to “on”. There is normally no need to adjust the voltage level. More information about this detector is available in the supplementary material.

The voltage output pulse of the detector is fed into a preamplifier attached to the PM tube, then to an amplifier in a NIM bin. [A Nuclear Instrument Module (NIM) bin is a rack with a number of slots for various instruments (amplifier, scaler, coincidence counter, etc) and a DC power supply for the modules. It allows great ease and flexibility in setting up an experiment.] After amplification, the signal is divided and fed into two separate pulse height analyzers, which measure the gamma ray spectrum in different ways.

Amplifier: The CANBERRA Model 2015A combines, in one single width module, a spectroscopy amplifier with gated baseline restoration and a timing single channel (SCA) pulse height analyzer. You will first use this unit as part of a multi-channel analyzer (MCA) and then as part of a single channel analyzer (SCA) and a coincidence circuit.

MCA function. In this part of the experiment the AMP output of the amplifier is fed into the MCA plug-in card and associated software on the computer. The MCA software sorts each PM pulse as it arrives from the amplifier according to pulse height and generates a fine-grained histogram of the number of voltage pulses from the detector versus pulse amplitude and, therefore, energy deposited in the scintillator. The amplifier output is proportional to the energy of the incoming signal with a maximum of 10V. The proportionality factor is controlled by the GAIN.

Note: If the gain is set too high the output corresponding to a signal of interest may be outside of the 10V range and the output will be saturated. The calibration procedure described below should prevent this.

SCA function. The SCA counts the number of events within a predetermined window of energies. In this mode the amplifier output is fed internally to the SCA, which has two knobs to adjust the range of pulse heights the SCA will pass.

The first knob adjusts the lower level E and the second pulse the window width ΔE , so that pulses with heights between E and $E + \Delta E$ will be passed to the output. The correspondence between energy and scale on the knobs is determined by the Gain, which you will set during the calibration. The SCA output is attached to the scaler input.

Scaler. The scaler counts the number of pulses arriving during a set time interval. To obtain the same data with the SCA as you would get with the MCA, you will need to step E across the range of possible voltages keeping ΔE constant and counting for a fixed length of time.

1. Calibration:

- a. Check that the apparatus is connected according to the diagram in Fig. 3 for detector 1. Turn on the NIM bin power supply and switch the high voltage PM supply to "standby". Turn on the computer. Log in as Student using the course password. The MCA program is called MAESTRO. Maximize the screen.
- b. Place the Na^{22} source on the Aluminum rod in front of detector 1. Turn the PM power supply from "standby" to "on" and make sure that the setting is about 1000 V (this would have been set by previous groups). Click on the "Go" button (in the top bar) to start accumulating the spectrum. It will take about five minutes to accumulate enough counts to resolve the spectrum. Note that the 1.2746 MeV photopeak is the highest energy feature in the spectrum, but not the strongest. Use the mouse to move the cursor onto the photopeak to identify its channel number (called "Marker" on the status bar at the bottom). Record the channel for the maximum of the photopeak. If it is in bin 900 skip the next step.
- c. Adjust the amplifier gain. Use the mouse to bring the cursor over the photopeak. The channel number and the number of counts will be listed at the bottom left of the screen. If the peak is not on channel 1900 adjust the gain accordingly and retake the spectrum. If necessary repeat a few times until the peak region falls in channel 1900 (plus/minus ~ 3). **Once you have set this gain do not change it;** otherwise you will change the calibration.
- d. To easily identify the energies of other features in the spectrum, click on the photopeak and go to "Calculate" \rightarrow "Calibration" and enter the energy of the photopeak in keV. Now, you can click on any other part of the spectrum and see the energy on the bottom of the screen.
- e. **Saving the data.** Go to "File" \rightarrow "Save As" and save your data in the ASCII SPE format. If you open this file in notepad, you will see a list of the counts in each bin. You can copy this into a data analysis program such as Origin or Excel.
- f. Now take the spectrum of the Cs^{137} source. Note the channel number of the 0.6616 MeV photopeak. Is it consistent with your calibration?
- g. Remove all sources to distant storage, and measure the background spectrum with the MCA.

2. **MCA Decay Spectra.** Without changing the calibrated gain setting record the spectra of the other two isotopes Cs^{137} , Co^{60} and Na^{22} with the MCA. *When taking data, be sure that you have removed the other sources far away from the scintillator!* When saving the data make sure that you include the collection time (you will need this later). For your report you will need to determine the (minimum) activity of each source in curies. To do this you have to measure (for each source) the appropriate distances to calculate the solid angle subtended by the scintillator relative to the source.

3. **SCA Decay Spectra.** Now measure the gamma spectrum of the cesium source using the SCA. Set the SCA window width at 0.1 V. [Note that there is a switch on the back of the module to change the full-scale for ΔE from 1 V to 10V.] Leave this setting fixed and step E from zero to a value beyond the Cs^{137} photopeak, recording at each step the scaler readings for a fixed count time (e.g. 100 s). In the region of the photopeak, where the counting rate changes rapidly with voltage, use 0.1 V steps in E. In the more slowly varying regions steps of 0.3 V should be adequate. Be sure to record all instrument settings, -- count time, PM voltage, etc. Plot the pulse height distribution spectrum [counting rate versus $E + (\Delta E/2)$].

4. **γ - γ correlation:** You will now study the correlation of the gamma rays emitted when the positron from the Na^{22} source annihilate. You will need two scintillation counters, two single channel analyzers (SCA), a coincidence module, and a scaler to observe this correlation. The first step will be to set the windows of the two SCAs to pass PM pulse heights corresponding to the 511 keV annihilation gamma rays. When a pulse from the PM tube meets the window requirements, the SCA outputs a narrow pulse, which is fed into the coincidence module. The coincidence module outputs a pulse whenever it detects pulses from both SCAs at the same time. The scaler counts these coincidence pulses for a fixed length of time. The two scintillators are mounted on arms that pivot about the location of the Na^{22} sample. When the two arms are 180° apart, strong coincidences will be recorded, but when the angle is decreased the count rate will drop sharply. (There will still be a number of accidental coincidences for nuclei that just happen to decay at the same time.)

4A. The first step will be to set up the window for counter 1. With the sodium source in place, observe the amplifier output with the MCA and with the oscilloscope (channel 1). [On the MCA spectrum of Na^{22} the 511 keV gamma ray is the strongest feature. Note that the positron annihilation occurs in the sodium source not in the detectors.] Set the scope sweep speed at $5\mu\text{s}/\text{cm}$ and the gain to 1 V/cm with the trigger source channel 1 and the trigger level about 2 V. You will see a series of superimposed pulses of widely varying height, but with many falling in the same band, which results in the scope traces being brightest for that

band. These pulses correspond to the 511 keV gamma rays. Adjust the amplifier gain so that this band corresponds to a pulse height of about 3 V.

Now connect the SCA output to channel 2, change the trigger source to channel 2, and set the scope to display both channels 1 and 2. [Note that the SCA pulses displayed in channel 2 are extremely narrow ($\sim 0.7 \mu\text{s}$.) which places stringent limits on whether the coincidence module will consider two gamma rays to be in coincidence.] Observe how the channel 1 scope trace changes as you change E, the lower edge of the window, and ΔE , the window width. Set E to about 2.8 V and ΔE to about 0.4V. [Note that there is a switch on the back of the module that allows you changes full-scale for ΔE from 1 V to 10V.] These settings will limit the pulses accepted by the SCA to just those corresponding to the 511 keV peak.

4B. Repeat step A for the second counter. Be sure to shift the cable leading to the MCA from the first amplifier to the second one.

4C. Now connect the outputs of the two SCAs to channels 1 and 2 of the coincidence module, set the corresponding switches to “coincidence” and connect the output of the coincidence module to the scaler input. Set the scaler to count for 10 or 100 s. With the two counters set 180° apart, you should obtain a strong coincidence count rate, which should drop almost to zero when the angle is reduced.

4D. Record the number of coincidences (for a fixed counting interval) as a function of the angle of separation of the detectors. Make a plot for your report. Explain the width of this curve, i.e. why do you obtain coincidences even when the angle is not exactly 180° ?

4E. For your report you will need to measure the appropriate distances to calculate the solid angle subtended by both scintillators relative to the source.

5. Co^{60} γ - γ correlation (optional): Figure 2 shows that when Co^{60} decays it emits two gamma rays. The second gamma ray rapidly follows the first (the excited state has a mean life of $0.7 \times 10^{-12}\text{s}$) and thus appears to be coincident. Its direction of emission is correlated with the direction in which the first is emitted. We will not go into the theory, which is complicated, but the prediction is that the angular correlation is given by

$$w(\theta) = 1 + \frac{\cos^2 \theta}{8} + \frac{\cos^4 \theta}{24}, \quad (2)$$

where θ is the angle between the two gamma rays. This is a rather small effect. The correlation function is 1.00 at 90° and 1.17 at 180° . This is a more challenging

measurement: to detect a 17% effect, you will need to measure the correlation rate to about 1%. The error in a measurement of N counts is \sqrt{N} . Thus to achieve 1% accuracy at a particular value of θ you will need to count until you have accumulated about 10^4 counts. For the source available in the lab this means you will need to count for 1000 s for each point.

5A. Mount the Co^{60} source at the center of the spectrometer and move the two detectors fairly close to the source (to get a large count rate).

5B. Set the windows for the two detectors to pass both of the gamma rays (1.1732 and 1.3325 MeV).

5C. Record the number of coincidences $C(\theta)$ as you vary θ in increments of 10° between 180° and 90° . Count for 1000 s for each point. Plot the ratio $C(\theta)/C(90^\circ)$ and compare with the theoretical correlation, Eq. (2).

ANALYSIS

Save each spectrum in a file and use your preferred software to analyze it. For each peak subtract the background count and record the net peak count, counting time, peak energy, FWHM and FWTM (full width at one-tenth of maximum) widths. On each plot label the photopeak(s) the Compton regions and all other features of interest. Be sure to study the k X-ray for the cesium source. Calculate the activity of each source in Curies. Make a plot of the

REPORT:

The report should contain the calibrated spectra for the sodium, cobalt, and cesium sources with labels identifying each spectral feature. It should contain a plot of the number of coincidences (for a fixed counting interval) as a function of the angle of separation of the detectors

1. From your data for cesium calculate the minimum activity in curies. Discuss why this can only be an estimate of the minimum activity. Make a spectral plot (counts versus $E + \Delta E/2$) of your SCA data for cesium and compare with the MCA spectrum.

2. For the cobalt spectrum determine the energies of the two gamma rays and compare with the accepted value. Discuss how the overlap of the tails affects the energies you measure. Calculate the theoretical Compton edges for the two gamma rays and indicate where they fall on your measured spectrum.

3. Resolution: From your MCA spectra calculate the resolution (FWHM) of the counter for the various photopeaks. Present your results in a table; is the

resolution dependent on the energy of the gamma ray?

4. Discuss all the spectral features you observe on the sodium spectrum. From the calibrated spectrum determine the rest mass of the electron with an estimate of your experimental error.
5. Make a plot of the correlation count rate versus angle between counters. Explain the width of the curve.
6. In the discussion session include answers to the following questions:
 - a. Why does the photoelectric peak not have zero width?
 - b. Why are there x-rays emitted after the ^{137}Cs decays? To what atomic transition and to what atom do these x-rays correspond?

REFERENCES:

experimental details

1. AN34: Experiments in Nuclear Science, 3rd Ed., EG & G Ortec. Experiment 1, Basic identifications in Electronic Measurement Systems, pp 1-7; Experiment 3, Gamma Spectroscopy Using NaI, pp 15-24; Linear and Logic Signals in EG & G Ortec NIM Instruments, pp159-160
2. D.W. Preston & E.R. Dietz: The Art of Experimental Physics, Appendix B Counting and Sorting Particles: The Scintillation Counter, pp376-385

general theory

3. Serway, Moses and Moyer: Modern Physics, pp 431-432 (scintillator/photomultiplier), 389-392 (beta and gamma decay)
4. Beiser: Concepts of Modern Physics, pp 443, 471-473 (gamma decay)
5. Halliday: Introductory Nuclear Physics
6. Siegbahn: Alpha, Beta and Gamma Ray Spectroscopy

JP 2/2000